

Product Information	
Catalog No.	PA11630
Product name	Anti-Brolucizumab Polyclonal Antibody
Applications	ELISA
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Rabbit polyclonal to Beovu.
Immunogen	Brolucizumab
Accession	CAS: 1531589-13-5
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	The product can be stored for 2 - 3 weeks at 2 to 8°C or for up to 12 months at -20°C. Avoid repeated freezing and thawing cycles.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Background

Brolucizumab is a humanized, single-chain variable fragment to the VEGF-A protein that blocks the effects of VEGF-A, a protein that promotes the growth of blood vessels in the choriocapillaris. Antibodysystem Anti-Brolucizumab Antibody, pAb, Rabbit is produced from sera of a rabbit immunized with Brolucizumab.