

Product Information	
Catalog No.	PA11636
Product name	Anti-Human CSNK2A1/CK2A1 Polyclonal Antibody
Applications	ELISA, WB
Species reactivity	Mouse, Bos taurus, Drosophila melanogaster, Chicken, Human, Oryctolagus cuniculus, Rat, Spodoptera frugiperda, Xenopus laevis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CSNK2A1, CK II alpha, CK2A1, Casein kinase II subunit alpha
Immunogen	E. coli - derived recombinant Human CSNK2A1/CK2A1 (Asp14-Ser337).
Accession	P68400, Q60737, Q8NEV1, P19139, P19784, O54833
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.