


Product Information	
Catalog No.	PA11712
Product name	Anti-Human TOMM20 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Danio rerio, Pongo abelii, Xenopus laevis, Xenopus tropicalis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Mitochondrial 20 kDa outer membrane protein, Mitochondrial import receptor subunit TOM20 homolog, KIAA0016, TOMM20, Outer mitochondrial membrane receptor Tom20
Immunogen	E. coli - derived recombinant Human TOMM20 (Arg51-Glu145).
Accession	Q15388
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-TOMM20 polyclonal antibody at 1µg/ml on various samples.

Lane 1: Hela cell lysate
Lane 2: HepG2 cell lysate
Lane 3: Jurkat cell lysate
Lane 4: THP-1 cell lysate