

Product Information	
Catalog No.	PA11767
Product name	Anti-Human UQCRC2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Bos taurus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	UQCRC2, Cytochrome b-c1 complex subunit 2, mitochondrial, Complex III subunit 2, Core protein II, Ubiquinol-cytochrome-c reductase complex core protein 2
Immunogen	E. coli - derived recombinant Human UQCRC2 (Pro35-Leu453).
Accession	P22695
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.