


Product Information	
Catalog No.	PA11833
Product name	Anti-Human SSBP1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Oryctolagus cuniculus, Pongo abelii
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Single-stranded DNA-binding protein, mitochondrial, Mt-SSB, MtSSB, PWP1-interacting protein 17, SSBP1, SSBP
Immunogen	E. coli - derived recombinant Human SSBP1 (His16-Glu148).
Accession	Q04837
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-SSBP1 polyclonal antibody at 1ug/ml on various samples.
Lane 1: A549 cell lysate
Lane 2: SW620 cell lysate