


Product Information	
Catalog No.	PA11842
Product name	Anti-Human TPI1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Dog, Chicken, Gorilla gorilla gorilla, Macaca fascicularis, Macaca mulatta, Mesocricetus auratus, Oryctolagus cuniculus, Pan troglodytes, Pongo abelii, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	TIM, Triose-phosphate isomerase, TPI, Methylglyoxal synthase, Triosephosphate isomerase, TPI1
Immunogen	E. coli - derived recombinant Human TPI1 (Met1-Gln249).
Accession	P60174
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-TPI1 polyclonal antibody at 1ug/ml on various samples.

- Lane 1: HepG2 cell lysate
- Lane 2: NIH3T3 cell lysate
- Lane 3: 293T cell lysate
- Lane 4: U87 cell lysate
- Lane 5: Jurkat cell lysate
- Lane 6: K562 cell lysate
- Lane 7: Hela cell lysate
- Lane 8: Mouse brain lysate
- Lane 9: Mouse liver lysate
- Lane 10: Mouse kidney lysate
- Lane 11: Mouse skeletal muscle lysate
- Lane 12: Rat brain lysate