

Product Information	
Catalog No.	PA11847
Product name	Anti-Human TNFAIP8L2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Callithrix jacchus, Oryctolagus cuniculus, Otolemur garnettii, Papio anubis, Rhinolophus ferrumequinum
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	TNF alpha-induced protein 8-like protein 2, Tumor necrosis factor alpha-induced protein 8-like protein 2, TNFAIP8-like protein 2, Inflammation factor protein 20, TIPE2, TNFAIP8L2
Immunogen	E. coli - derived recombinant Human TNFAIP8L2 (Met1-Leu184).
Accession	Q6P589
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.