

Product Information	
Catalog No.	PA11889
Product name	Anti-Human MDA5/IFIH1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	MDA-5, IFIH1, Helicard, Clinically amyopathic dermatomyositis autoantigen 140 kDa, Interferon-induced with helicase C domain protein 1, Helicase with 2 CARD domains, RNA helicase-DEAD box protein 116, RIG-I-like receptor 2, Interferon-induced helicase C domain-containing protein 1, MDA5, Murabutide down-regulated protein, Melanoma differentiation-associated protein 5, CADM-140 autoantigen, RH116, RLR-2
Immunogen	E. coli - derived recombinant Human MDA5/IFIH1 (Arg309-Asp1025).
Accession	Q9BYX4
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.