

Product Information	
Catalog No.	PA11892
Product name	Anti-Human WDR1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Chicken, Xenopus laevis, Xenopus tropicalis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	WD repeat-containing protein 1, Actin-interacting protein 1, WDR1, AIP1, NORI-1
Immunogen	E. coli - derived recombinant Human WDR1 (Met1-Lys317).
Accession	O75083
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.