

Product Information	
Catalog No.	PA11897
Product name	Anti-Human MICAL1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse, Bos taurus, Human, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	NICAL, [F-actin]-monooxygenase MICAL1, MICAL1, MICAL-1, NEDD9-interacting protein with calponin homology and LIM domains, Molecule interacting with CasL protein 1, MICAL
Immunogen	E. coli - derived recombinant Human MICAL1 (Met1-Lys488).
Accession	Q8TDZ2
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.