

Product Information	
Catalog No.	PA11905
Product name	Anti-Human TCF4 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Dog, Rat, Danio rerio, Chicken, Mesocricetus auratus, Xenopus laevis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	TCF-4, SL3-3 enhancer factor 2, Immunoglobulin transcription factor 2, ITF2, SEF-2, SEF2, bHLHB19, ITF-2, TCF4, Transcription factor 4, BHLHB19, Class B basic helix-loop-helix protein 19
Immunogen	E. coli - derived recombinant Human TCF4 (Asp552-Ala631).
Accession	P15884
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.