

Product Information	
Catalog No.	PA11937
Product name	Anti-Human PRKCB Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Danio rerio, Oryctolagus cuniculus, Xenopus laevis, Xenopus tropicalis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	PKC-beta, PRKCB1, PKC-B, PKCB, PRKCB, Protein kinase C beta type
Immunogen	E. coli - derived recombinant Human PRKCB (Asp341-Asn664).
Accession	P05771
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.