

Product Information	
Catalog No.	PA11965
Product name	Anti-Human CDC6 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	HsCdc18, CDC18L, CDC6, HsCDC6, CDC6-related protein, Cell division control protein 6 homolog, p62(cdc6), Cdc18-related protein
Immunogen	E. coli - derived recombinant Human CDC6 (Lys151-Pro560).
Accession	Q99741
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.