

Product Information	
Catalog No.	PA11970
Product name	Anti-Human TRIM28 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	TIF1-beta, Tripartite motif-containing protein 28, KRAB-associated protein 1, TRIM28, RING-type E3 ubiquitin transferase TIF1-beta, KRAB-interacting protein 1, KRIP-1, E3 SUMO-protein ligase TRIM28, RNF96, Nuclear corepressor KAP-1, TIF1B, Transcription intermediary factor 1-beta, KAP1, KAP-1, RING finger protein 96
Immunogen	E. coli - derived recombinant Human TRIM28 (Ser50-Arg136).
Accession	Q13263
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.