

Product Information	
Catalog No.	PA12002
Product name	Anti-Human ATF6B Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CREBL1, Protein G13, cAMP-responsive element-binding protein-like 1, Activating transcription factor 6 beta, cAMP response element-binding protein-related protein, Cyclic AMP-dependent transcription factor ATF-6 beta, G13, ATF6B, Creb-rp, ATF6-beta, cAMP-dependent transcription factor ATF-6 beta
Immunogen	E. coli - derived recombinant Human ATF6B (Ala2-Arg50).
Accession	Q99941
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.