

Product Information	
Catalog No.	PA12020
Product name	Anti-Human PINX1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Pin2-interacting protein X1, LPTS, TRF1-interacting protein 1, Liver-related putative tumor suppressor, PIN2/TERF1-interacting telomerase inhibitor 1, LPTL, PINX1, Protein 67-11-3
Immunogen	E. coli - derived recombinant Human PINX1 (Asn72-Lys149).
Accession	Q96BK5
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.