

Product Information	
Catalog No.	PA12021
Product name	Anti-Human PRKACA Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Dog, Rat, Bos taurus, Caenorhabditis briggsae, Caenorhabditis elegans, Cricetulus griseus, Drosophila melanogaster, Macaca mulatta, Ovis aries, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	PRKACA, PKACA, PKA C-alpha, cAMP-dependent protein kinase catalytic subunit alpha
Immunogen	E. coli - derived recombinant Human PRKACA (His40-Lys310).
Accession	P17612
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.