


Product Information	
Catalog No.	PA12028
Product name	Anti-Human RSAD2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus, Pongo abelii, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	RSAD2, Viperin, Cytomegalovirus-induced gene 5 protein, Radical S-adenosyl methionine domain-containing protein 2, Virus inhibitory protein, endoplasmic reticulum-associated, interferon-inducible, CIG5
Immunogen	E. coli - derived recombinant Human RSAD2 (Ser74-Trp361).
Accession	Q8WYG1
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-RSAD2 polyclonal antibody at 1ug/ml on various samples.

- Lane 1: HUVEC cell lysate
- Lane 2: Mouse spleen lysate
- Lane 3: Mouse brain lysate
- Lane 4: Mouse lung lysate
- Lane 5: Rat spleen lysate
- Lane 6: Rat brain lysate
- Lane 7: Rat lung lysate