


Product Information	
Catalog No.	PA12052
Product name	Anti-Human IPO4 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	RanBP4, Imp4b, Importin-4b, Ran-binding protein 4, RANBP4, IPO4, Importin-4, Imp4, IMP4B
Immunogen	E. coli - derived recombinant Human IPO4 (Asn667-Ser1081).
Accession	Q8TEX9
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-IPO4 polyclonal antibody at 1ug/ml on various samples.

Lane 1: Hela cell lysate

Lane 2: Lncap cell lysate

Lane 3: A549 cell lysate