


Product Information	
Catalog No.	PA12087
Product name	Anti-Human ATP1B1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Sodium/potassium-transporting ATPase subunit beta-1, ATP1B1, Sodium/potassium-dependent ATPase subunit beta-1, ATP1B
Immunogen	E. coli - derived recombinant Human ATP1B1 (Gly76-Ser303).
Accession	P05026
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-ATP1B1 polyclonal antibody at 1ug/ml on various samples.
Lane 1: Mouse brain lysate
Lane 2: Rat brain lysate