

Product Information	
Catalog No.	PA12138
Product name	Anti-Citrus clementina Laccase Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Citrus clementina (Clementine) (Citrus deliciosa x Citrus sinensis)
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Laccase, 1.10.3.2, Benzenediol:oxygen oxidoreductase, Diphenol oxidase, Urishiol oxidase, CICLE_v10031134mg
Immunogen	E. coli - derived recombinant Citrus clementina Laccase (Ala22-Cys553).
Accession	V4SM97
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.