


Product Information	
Catalog No.	PA12159
Product name	Anti-Human FEN1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	FEN-1, Maturation factor 1, hFEN-1, DNase IV, RAD2, FEN1, Flap endonuclease 1, MF1, Flap structure-specific endonuclease 1
Immunogen	E. coli - derived recombinant Human FEN1 (Met1-Thr336).
Accession	P39748
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-FEN1 polyclonal antibody at 1ug/ml on various samples.

- Lane 1: Daudi cell lysate
- Lane 2: NIH3T3 cell lysate
- Lane 3: HepG2 cell lysate
- Lane 4: 293T cell lysate
- Lane 5: Jurkat cell lysate
- Lane 6: PC-12 cell lysate
- Lane 7: Mouse testis lysate
- Lane 8: Mouse spleen lysate
- Lane 9: Rat testis lysate
- Lane 10: Rat spleen lysate