

Product Information	
Catalog No.	PA12273
Product name	Anti-Human CELF2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CELF-2, hNAPOR, Neuroblastoma apoptosis-related RNA-binding protein, CUGBP Elav-like family member 2, CUG-BP- and ETR-3-like factor 2, BRUNOL3, NAPOR, ETR3, CUGBP2, CUG triplet repeat RNA-binding protein 2, ETR-3, CELF2, ELAV-type RNA-binding protein 3, RNA-binding protein BRUNOL-3, Bruno-like protein 3, CUG-BP2
Immunogen	E. coli - derived recombinant Human CELF2 (Lys41-Asp211).
Accession	O95319
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.