

Product Information	
Catalog No.	PA12304
Product name	Anti-Human RBCK1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	RING finger protein 54, RNF54, RanBP-type and C3HC4-type zinc finger-containing protein 1, C20orf18, XAP4 X-associated protein 4, UBCE7IP3, RBCK1, RING-type E3 ubiquitin transferase HOIL-1, HOIL-1, Heme-oxidized IRP2 ubiquitin ligase 1, XAP3, Ubiquitin-conjugating enzyme 7-interacting protein 3, HBV-associated factor 4
Immunogen	E. coli - derived recombinant Human RBCK1 (Met1-Gly156).
Accession	Q9BYM8
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.