


Product Information	
Catalog No.	PA12360
Product name	Anti-Human TXNL1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Thioredoxin-like protein 1, TXNL, TXNL1, 32 kDa thioredoxin-related protein, TXL, TRP32
Immunogen	E. coli - derived recombinant Human TXNL1 (Met1-Gln271).
Accession	O43396
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-TXNL1 polyclonal antibody at 1ug/ml on various samples.

- Lane 1: K562 cell lysate
- Lane 2: HL-60 cell lysate
- Lane 3: SKOV3 cell lysate
- Lane 4: Caco2 cell lysate
- Lane 5: NIH3T3 cell lysate
- Lane 6: PC-12 cell lysate
- Lane 7: HCT-116 cell lysate
- Lane 8: Mouse brain lysate
- Lane 9: Mouse heart lysate
- Lane 10: Rat brain lysate
- Lane 11: Rat heart lysate