

Product Information	
Catalog No.	PA12403
Product name	Anti-Human ZFP36L1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	ERF-1, ZFP36-like 1, TIS11B, BRF1, ZFP36L1, Butyrate response factor 1, Zinc finger protein 36, C3H1 type-like 1, ERF1, RNF162B, BERG36, TPA-induced sequence 11b, mRNA decay activator protein ZFP36L1, EGF-response factor 1
Immunogen	E. coli - derived recombinant Human ZFP36L1 (Arg121-Arg195).
Accession	Q07352
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.