

Product Information	
Catalog No.	PA12404
Product name	Anti-Human ZKSCAN3 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	ZNF309, Zinc finger protein 47 homolog, ZKSCAN3, Zfp-47, Zinc finger protein 306, ZSCAN13, ZNF306, Zinc finger protein 309, Zf47, ZFP47, Zinc finger protein with KRAB and SCAN domains 3, Zinc finger and SCAN domain-containing protein 13
Immunogen	E. coli - derived recombinant Human ZKSCAN3 (Gly41-Gly176).
Accession	Q9BRR0
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.