

Product Information	
Catalog No.	PA12442
Product name	Anti-Human BHLHE40 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	DEC1, bHLHe40, Differentially expressed in chondrocytes protein 1, Class E basic helix-loop-helix protein 40, Class B basic helix-loop-helix protein 2, STRA13, SHARP2, SHARP-2, bHLHb2, Enhancer-of-split and hairy-related protein 2, BHLHB2, BHLHE40, Stimulated by retinoic acid gene 13 protein
Immunogen	E. coli - derived recombinant Human BHLHE40 (Ser46-Gly186).
Accession	O14503
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.