

Product Information	
Catalog No.	PA12503
Product name	Anti-Human PIP5K1A Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	PIP5K1-alpha, PtdIns(4)P-5-kinase 1 alpha, Phosphatidylinositol 4-phosphate 5-kinase type-1 alpha, PIP5K1alpha, 68 kDa type I phosphatidylinositol 4-phosphate 5-kinase alpha, PIP5K1A, Phosphatidylinositol 4-phosphate 5-kinase type I alpha
Immunogen	E. coli - derived recombinant Human PIP5K1A (Lys62-Lys452).
Accession	Q99755
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.