

Product Information	
Catalog No.	PA12597
Product name	Anti-Human SSRP1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Recombination signal sequence recognition protein 1, SSRP1, FACT 80 kDa subunit, Chromatin-specific transcription elongation factor 80 kDa subunit, T160, FACT complex subunit SSRP1, FACTp80, FACT80, Structure-specific recognition protein 1, hSSRP1, Facilitates chromatin transcription complex subunit SSRP1, Facilitates chromatin transcription complex 80 kDa subunit
Immunogen	E. coli - derived recombinant Human SSRP1 (Met1-Lys422).
Accession	Q08945
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.