

Product Information	
Catalog No.	PA12615
Product name	Anti-Human VPS4B Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	SKD1, VPS4B, Suppressor of K(+) transport growth defect 1, Cell migration-inducing gene 1 protein, Protein SKD1, Vacuolar protein sorting-associated protein 4B, VPS42
Immunogen	E. coli - derived recombinant Human VPS4B (Gly138-Gly444).
Accession	O75351
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.