

Product Information	
Catalog No.	PA12662
Product name	Anti-Human TRIM9 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Tripartite motif-containing protein 9, RING-type E3 ubiquitin transferase TRIM9, RNF91, RING finger protein 91, TRIM9, KIAA0282, E3 ubiquitin-protein ligase TRIM9
Immunogen	E. coli - derived recombinant Human TRIM9 (Asp541-Ala710).
Accession	Q9C026
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.