

Product Information	
Catalog No.	PA12688
Product name	Anti-Human MEPCE Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Dog, Cat, Bovine, Pig, Alpaca, Camel, Ferret, Horse, Lion, Rabbit, Red fox, Sheep, Tree shrew
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	MEPCE, Bin3 homolog, Bicoid-interacting protein 3 homolog, 7SK snRNA methylphosphate capping enzyme, MePCE, BCDIN3
Immunogen	E. coli - derived recombinant Human MEPCE (Lys407-His689).
Accession	Q7L2J0
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.