

Product Information	
Catalog No.	PA12718
Product name	Anti-Human ACTL6A Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	BAF53, BAF53A, INO80K, BRG1-associated factor 53A, ArpNbeta, Actin-related protein Baf53a, INO80 complex subunit K, ACTL6A, Actin-like protein 6A, 53 kDa BRG1-associated factor A
Immunogen	E. coli - derived recombinant Human ACTL6A (Gly12-Pro429).
Accession	O96019
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.