

Product Information	
Catalog No.	PA12722
Product name	Anti-Human IFNGR2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	IFN-gamma receptor 2, Interferon gamma receptor beta-chain, Interferon gamma transducer 1, IFN-gamma-R-beta, AF-1, IFN-gamma-R2, IFNGT1, Interferon gamma receptor accessory factor 1, Interferon gamma receptor 2, IFNGR2
Immunogen	E. coli - derived recombinant Human IFNGR2 (Lys37-Asp240).
Accession	P38484
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.