

Product Information	
Catalog No.	PA12734
Product name	Anti-Human CEP290 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Nephrocystin-6, KIAA0373, CEP290, Bardet-Biedl syndrome 14 protein, Tumor antigen se2-2, BBS14, NPHP6, Centrosomal protein of 290 kDa, Cancer/testis antigen 87, Cep290, CT87
Immunogen	E. coli - derived recombinant Human CEP290 (Thr1289-Arg1416).
Accession	O15078
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.