

Product Information	
Catalog No.	PA12736
Product name	Anti-Human MRPL58 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	DS1, MRPL58, Immature colon carcinoma transcript 1 protein, Digestion substraction 1, Mitochondrial large ribosomal subunit protein mL62, DS-1, ICT1, MRP-L58, Peptidyl-tRNA hydrolase ICT1, mitochondrial, 39S ribosomal protein L58, mitochondrial
Immunogen	E. coli - derived recombinant Human MRPL58 (Asp71-Asp206).
Accession	Q14197
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.