

Product Information	
Catalog No.	PA12821
Product name	Anti-BAIAP2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	BAI-associated protein 2, Insulin receptor substrate p53, IRS-58, BAI1-associated protein 2, IRSp53, Insulin receptor substrate p53/p58, Fas ligand-associated factor 3, BAIAP2, FLAF3, Insulin receptor substrate protein of 53 kDa, Brain-specific angiogenesis inhibitor 1-associated protein 2, Protein BAP2, IRSp53/58
Immunogen	E. coli - derived recombinant Human BAIAP2 (Met1-Asn300).
Accession	Q9UQB8
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.