


Product Information	
Catalog No.	PA12849
Product name	Anti-SP1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	SP1, TSFP1, Transcription factor Sp1
Immunogen	E. coli - derived recombinant Human SP1 (Leu260-Ala591).
Accession	P08047
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-SP1 polyclonal antibody at 1ug/ml on various samples.

Lane 1: Jurkat cell lysate

Lane 2: K562 cell lysate

Lane 3: HepG2 cell lysate