

Product Information	
Catalog No.	PA12987
Product name	Anti-ABCD4 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	ATP-binding cassette sub-family D member 4, ABCD4, PXMP1-L, PMP70-related protein, Peroxisomal membrane protein 1-like, P70R, PMP69, Peroxisomal membrane protein 69, PXMP1L, Lysosomal cobalamin transporter ABCD4
Immunogen	E. coli - derived recombinant Human ABCD4 (Leu356-Glu606).
Accession	O14678
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.