

Product Information	
Catalog No.	PA12989
Product name	Anti-FBLN5 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	UP50, Urine p50 protein, DANCE, Fibulin-5, Developmental arteries and neural crest EGF-like protein, Dance, FBLN5, FIBL-5
Immunogen	E. coli - derived recombinant Human FBLN5 (Leu99-Cys205).
Accession	Q9UBX5
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.