

Product Information	
Catalog No.	PA13035
Product name	Anti-KDM4A Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Lysine-specific demethylase 4A, JMJD2, Jumonji domain-containing protein 2A, KDM4A, KIAA0677, JHDM3A, [histone H3]-trimethyl-L-lysine(9) demethylase 4A, [histone H3]-trimethyl-L-lysine(36) demethylase 4A, JMJD2A, JmjC domain-containing histone demethylation protein 3A
Immunogen	E. coli - derived recombinant Human KDM4A (Tyr141-Lys310).
Accession	O75164
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.