

Product Information	
Catalog No.	PA13128
Product name	Anti-KCNJ10 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	ATP-dependent inwardly rectifying potassium channel Kir4.1, KCNJ10, ATP-sensitive inward rectifier potassium channel 10, Inward rectifier K(+) channel Kir1.2, Potassium channel, inwardly rectifying subfamily J member 10
Immunogen	E. coli - derived recombinant Human KCNJ10 (Phe165-Val379).
Accession	P78508
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.