

Product Information	
Catalog No.	PA13178
Product name	Anti-CPT1B Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CPT1B, Carnitine palmitoyltransferase 1B, CPT1-M, Carnitine O-palmitoyltransferase I, muscle isoform, CPT I, Carnitine O-palmitoyltransferase 1, muscle isoform, CPT1-M, Carnitine palmitoyltransferase I-like protein, KIAA1670
Immunogen	E. coli - derived recombinant Human CPT1B (Thr602-Ser772).
Accession	Q92523
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.