


Product Information	
Catalog No.	PA13179
Product name	Anti-ANXA7 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	ANX7, Annexin A7, SNX, Annexin VII, Synexin, Annexin-7, ANXA7
Immunogen	E. coli - derived recombinant Human ANXA7 (Met1-Lys287).
Accession	P20073
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-ANXA7 polyclonal antibody at 1ug/mL on recombinant ANXA7 (Catalog No. HB127012).

Western blot