

Product Information	
Catalog No.	PA13272
Product name	Anti-Human ATP5IF1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	IF1, Inhibitor of F(1)F(o)-ATPase, ATP synthase F1 subunit epsilon, IF(1), ATP5IF1, ATPi, ATPiF1, ATPase inhibitor, mitochondrial
Immunogen	E. coli - derived recombinant Human ATP5IF1 (Met1-Asp106).
Accession	Q9UII2
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.