

Product Information	
Catalog No.	PA13305
Product name	Anti-CD315/PTGFRN Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CD9 partner 1, PTGFRN, EWI-F, CD9P-1, CD315, CD9P1, FPRP, Prostaglandin F2 receptor negative regulator, Prostaglandin F2-alpha receptor-associated protein, KIAA1436, Glu-Trp-Ile EWI motif-containing protein F, Prostaglandin F2-alpha receptor regulatory protein, EWIF
Immunogen	E. coli - derived recombinant Human CD315/PTGFRN (Pro406-Pro662).
Accession	Q9P2B2
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.