

Product Information	
Catalog No.	PA13309
Product name	Anti-TNKS2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	TANK2, ADP-ribosyltransferase diphtheria toxin-like 6, Tankyrase-like protein, Poly [ADP-ribose] polymerase tankyrase-2, Tankyrase-related protein, TNKS2, Tankyrase II, ARTD6, Poly [ADP-ribose] polymerase 5B, Tankyrase-2, PARP5B, TNKL, TRF1-interacting ankyrin-related ADP-ribose polymerase 2, TNKS-2, Protein poly-ADP-ribosyltransferase tankyrase-2
Immunogen	E. coli - derived recombinant Human TNKS2 (Ser959-Gly1166).
Accession	Q9H2K2
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.