


Product Information	
Catalog No.	PA13351
Product name	Anti-CD56/NCAM1 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	NCAM-1, NCAM1, CD56, N-CAM-1, NCAM, Neural cell adhesion molecule 1
Immunogen	E. coli - derived recombinant Human CD56/NCAM1 (Glu30-Glu352).
Accession	P13591
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-NCAM1 polyclonal antibody at 1ug/ml on various samples.

Lane 1: Mouse cerebellum lysate

Lane 2: Rat brain lysate

Lane 3: Rat cerebellum lysate