

Product Information	
Catalog No.	PA13353
Product name	Anti-MAPK12 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Mitogen-activated protein kinase 12, ERK6, MAPK12, Extracellular signal-regulated kinase 6, ERK-6, MAP kinase p38 gamma, SAPK3, MAPK 12, Stress-activated protein kinase 3, MAP kinase 12, Mitogen-activated protein kinase p38 gamma
Immunogen	E. coli - derived recombinant Human MAPK12 (Tyr27-Phe311).
Accession	P53778
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.